

Disposable Cytology Brush

Remark: Special specification can be customized.

Disposable Cytology Brush

Straight type, with sheath.



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRS-Aφ 1.8*2300-PN	φ1.8x2300	≥2.0	VET-6010, VET-8015, MVE-6010, MVE-5010

Cleaning Brush

Remark: Special specification can be customized.

Disposable Cleaning Brush



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRX1823	φ1.8x2300	≥2.0	VET-6010, VET-8015, MVE-6010, MVE-5010
JRX2423	φ2.4x2300	≥2.8	VET-9215
JRX2440	φ2.4x4000	≥2.8	VET-9830, VET-1330, VET-1335

Reusable Cleaning Brush

Autoclavable 134 C



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRX1823	φ1.8x2300	≥2.0	VET-6010, VET-8015, MVE-6010, MVE-5010

Stent System



Magellan Medical Systems

CE 0197
ISO 9001/ ISO 13485 Certified

Web: www.magellanmedsys.com

Email: support@magellanmedsys.com

Tel: 949-878-1690

AOHUA

AOHUA



Veterinary
Endoscopic Accessories Solution

Disposable Biopsy Forceps

Collect tissue from human body through endoscope.

Disposable Oval Cup Biopsy Forceps

Oval, large opening cups for gathering bigger tissue sample.



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRQ-Y1823-HPA	φ1.8x2300	≥2.0	VET-6010, VET-8015, MVE-6010, MVE-5010
JRQ-Y2323-HPA	φ2.3x2300	≥2.8	VET-9215
JRQ-Y2340-HPA	φ2.3x4000	≥2.8	VET-9830, VET-1330, VET-1335



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRQ-Y1823-HPAC	φ1.8x2300	≥2.0	VET-6010, VET-8015, MVE-6010, MVE-5010
JRQ-Y2323-HPAC	φ2.3x2300	≥2.8	VET-9215
JRQ-Y2340-HPAC	φ2.3x4000	≥2.8	VET-9830, VET-1330, VET-1335

Disposable Oval Cup with Needle Biopsy Forceps

Oval, large opening cups for gathering bigger tissue sample. With needle for accurate targeting, no sliding.



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRQ-Z2323-HPA	φ2.3x2300	≥2.8	VET-9215
JRQ-Z2340-HPA	φ2.3x4000	≥2.8	VET-9830, VET-1330, VET-1335



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRQ-Z2323-HPAC	φ2.3x2300	≥2.8	VET-9215
JRQ-Z2340-HPAC	φ2.3x4000	≥2.8	VET-9830, VET-1330, VET-1335

Disposable Alligator Cup Biopsy Forceps

Alligator cups make the acquisition of smooth tissue.



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRQ-Y2323-HEA	φ2.3x2300	≥2.8	VET-9215
JRQ-Y2340-HEA	φ2.3x4000	≥2.8	VET-9830, VET-1330, VET-1335



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRQ-Y2321-HEAC	φ2.3x2300	≥2.8	VET-9215
JRQ-Y2340-HEAC	φ2.3x4000	≥2.8	VET-9830, VET-1330, VET-1335

Remark:

- ◆ Working Length, Diameter of Biopsy Forceps can be customized.
- ◆ Optional Plastic Coating Layer, which can effectively prevent the biopsy forceps hurting the Endoscope channel.

Disposable Retrieval Basket

Remark: Special specification can be customized.

Disposable Retrieval Basket



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRY-WD-1821-25	φ1.8x2100	≥2.0	VET-8015
JRY-WD-2421-35	φ2.4x2100	≥2.8	VET-9215
JRQ-WD-1815-25	φ1.8x1500	≥2.0	VET-6010, VME-6010, VME-5010

Disposable Snare

Remark: Special specification can be customized.

Disposable Snare



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRY-QA-1821-25	φ1.8x2100	≥2.0	VET-8015
JRY-QA-2421-30	φ2.4x2100	≥2.8	VET-9215
JRQ-QA-1815-25	φ1.8x1500	≥2.0	VET-6010, VME-6010, VME-5010

Disposable Grasping Forceps

Remark: Special specification can be customized.

Rat Tooth Forceps



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRY-QD-1821	φ1.8x2100	≥2.0	VET-8015
JRY-QD-2421	φ2.4x2100	≥2.8	VET-9215
JRY-QD-1815	φ1.8x1500	≥2.0	VET-6010, VME-6010, VME-5010

Pelican Type Grasping Forceps



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRY-QT-2421	φ2.4x2100	≥2.8	VET-9215

Alligator Grasper



Item No.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Compatible
JRY-QE-2421	φ2.4x2100	≥2.8	VET-9215